


FEHLING...

... THE DIFFERENCE

8


SUCTION INSTRUMENTS


SAUGINSTRUMENTE



FEHLING
INSTRUMENTS

Hanauer Landstraße 7A · 63791 Karlstein/Germany · www.fehling-instruments.de

 +49 (0) 61 88 - 95 74.40

 +49 (0) 61 88 - 95 74.45

 info@fehling-instruments.de



SUPERPLAST
A = 180 mm

Ø TUBE / ROHR	
2,0 mm	HAY-2
3,0 mm	HAY-3
4,0 mm	HAY-4

SUPERPLAST instruments are malleable at regular OR temperature. They can be deformed to any helpful shape by soft bending power. While exposed to higher temperatures as i.e. that one of a washing machine or an autoclave, they automatically regain their original shape.


SUPERPLAST Instrumente sind bei normaler OP-Temperatur biegsam. Sie können durch sanften Druck in die gewünschte Form gebracht werden. Unter höheren Temperaturen, wie bei der Sterillisation, nehmen sie ihre ursprüngliche Form wieder an.



MALLEABLE SUCTION CANNULAS
BIEGSAME SAUGER


Ø TUBE / ROHR	LENGTH / LÄNGE A						
	80 mm	100 mm	125 mm	145 mm	165 mm	180 mm	220 mm
0,6 mm	HAU-0	-	-	-	-	-	-
0,8 mm	HAU-4	-	-	HAI-0	-	-	-
1,0 mm	HAU-5	HAU-1	HAU-6	HAI-7	HAX-4	-	-
1,2 mm	-	HAU-2	HAU-7	HAI-8	-	-	-
1,6 mm	-	HAU-3	HAU-8	HAI-9	-	-	-
1,8 mm	-	-	-	-	HAI-1	HAI-4	HAX-1
2,0 mm	-	-	-	HAU-9	HAI-2	HAI-5	HAX-2
2,3 mm	-	-	-	-	HAI-3 (2,2 mm)	HAI-6	HAX-3
3,0 mm	-	-	-	-	HAX-5	-	-

FRAZIER




CERAMO®

Ø TUBE / ROHR	LENGTH / LÄNGE A / B		
	115 / 180 mm	145 / 200 mm	210 / 250 mm
1,0 mm	–	VOC-1	VOE-4
2,0 mm	VOB-2	VOC-2	VOE-5
2,5 mm	VOB-3	VOC-3	VOE-6
3,0 mm	VOB-4	VOC-4	VOE-7
4,0 mm	VOB-5	VOC-5	VOE-8
5,0 mm	VOB-6	VOC-6	VOE-9



UNCOATED
UNBESCHICHTET

Ø TUBE / ROHR	LENGTH / LÄNGE A / B		
	115 / 180 mm	145 / 200 mm	210 / 250 mm
1,0 mm	–	HAR-1	–
2,0 mm	NAT-7	HAR-2	NHR-8
2,5 mm	NAT-8	HAR-0	NHS-0
3,0 mm	NAT-9	HAR-3	NHR-9
4,0 mm	NAU-1	HAR-4	NHS-1
5,0 mm	NAU-2	HAP-9	NHS-2





CERAMO®

MALLEABLE / BIEGBAR

Ø TUBE / ROHR	LENGTH / LÄNGE A / B	
	120 / 180 mm	160 / 210 mm
1,0 mm	VOA-5	VOB-0
2,0 mm	VOA-6	VOA-1
3,0 mm	VOA-7	VOA-2
4,0 mm	VOA-8	VOA-3
5,0 mm	VOA-9	VOA-4
6,0 mm	-	VOA-0




Ø TIP/SPITZE	
 3,3 mm	WEU-0
 4,0 mm	WEU-1

NON-STICK
MONOPOLAR SUCTION TUBE
MONOPOLARER SAUGER
A = 130 mm


KINTZEL

DISSECTION CANNULAS
DISSEKTIONSSAUGER



CERAMO®

Ø TUBE / ROHR	LENGTH / LÄNGE A / B	
3,0 mm	125 / 190 mm	HAP-5
4,0 mm	175 / 245 mm	HAP-4



UNCOATED
UNBESCHICHTET

Ø TUBE / ROHR	LENGTH / LÄNGE A / B	
3,0 mm	110 / 170 mm	NHR-0
4,0 mm	120 / 190 mm	NHR-4
4,0 mm	140 / 210 mm	NHA-8

JAMIESON

DISSECTION CANNULAS
DISSEKTIONSSAUGER



	HOLE/BALL LOCH/KUGEL Ø	LENGTH LÄNGE	
	0,5 mm	270 mm	MMT-1
	1,0 mm	210 mm	MNF-0
	1,5 mm	270 mm	MMA-7
	1,0 mm	260 mm	MMC-0
	1,5 mm	290 mm	MME-1
	1,5 mm	290 mm	MME-2
	0,5 mm	270 mm	MMJ-0
	1,0 mm	270 mm	MMB-0
	1,0 mm	270 mm	MMJ-6



VALVULAR SUCTION TUBE
KLAPPENSAUGER
300 mm
MNO-0



CERAMO®
SUCTION CANNULA FOR
VALVE DECALCIFICATION
ENTKALKUNGSSAUGER
240 mm
MNF-5



HLM-ASPIRATOR
HLM-SAUGER
280 mm
MEL-7

PEDIATRIC
PÄDIATRISCH



225 mm
MCE-9



245 mm
MCE-8



210 mm
MNE-5



SEBENING
180 mm
MNE-8

GRECO

PEDIATRIC
PÄDIATRISCH



DETACHABLE STRAINER
KORB ABNEHMBAR

180 mm

∅ TUBE / ROHR	
3,3 mm	MNM-1
5,0 mm	MNM-3



DETACHABLE STRAINER
KORB ABNEHMBAR

160 mm

∅ TUBE / ROHR	
3,3 mm	MNM-2
5,0 mm	MNM-4



210 mm

∅ TUBE / ROHR	
3,3 mm	MNL-2
5,0 mm	MNL-5



180 mm

∅ TUBE / ROHR	
3,3 mm	MNL-1
5,0 mm	MNL-4





150 mm

∅ TUBE / ROHR	
3,3 mm	MNL-3
5,0 mm	MNL-6

GRECO


PEDIATRIC
PÄDIATRISCH



Ø TUBE / ROHR	
3,3 mm	MNM-5
5,0 mm	MNM-7

Ø TUBE / ROHR	
3,3 mm	MNM-6
5,0 mm	MNM-8

VENTRICULAR SUCTION CANNULAS
VENTRIKELSAUGER



Ø 6 mm
L 23 mm
MKN-8

Ø 7 mm
L 25 mm
MKK-5

1/4

PEDIATRIC
PÄDIATRISCH



FOR 1/4" TUBE
FÜR 1/4" SCHLAUCH


Ø 16 mm
MCC-0

1/4




9 x 12 mm
MCC-1

1/4





Ø 2,5 mm L 48 mm
HAH-9
Ø 3,7 mm L 55 mm
HAJ-2
Ø 5,0 mm L 160 mm
HAJ-3

HEPARIN CANNULA
HEPARINKANÜLE




MALLEABLE / BIEGSAM

Ø	1,2 mm	1,5 mm	2,0 mm
	VAD-3	VAD-5	VAD-0
	VAD-4	VAD-6	

BUTTONED CANNULAS
KNOPFKANÜLEN
80 mm

VASCULAR IRRIGATORS
VASKULARSPÜLKÖPFE



LUER-LOCK ADAPTER
LUER-LOCK-ANSATZ