


FEHLING...
... THE DIFFERENCE

16

ROTARY
INSTRUMENTS
ROTIERENDE
INSTRUMENTE



FEHLING
INSTRUMENTS

Hanauer Landstraße 7A · D-63791 Karlstein · www.fehling-instruments.de



+49 (0) 61 88 - 9574.0




+49 (0) 61 88 - 9574.44




info@fehling-instruments.de

Tables and Designations / Tabellen und Bezeichnungen


LOW SPEED Cutters / Fräser























HIGH SPEED Cutters / Fräser





























LOW SPEED Twist Drills / Spiralbohrer












CERAMO® LOW SPEED Rose Hobs / Rosenbohrer

							
	Ø mm	S	L 70 mm	L 95 mm	L 125 mm	L 55 mm	L 112 mm
	0,6	6	LTK-8	LRK-6	LTE-9		
	0,8	6	LRI-2	LTA-0	LRG-1		
		12	LTK-9				
	1,0	6	LTL-1	LTA-1	LRG-2		
		12	LRI-4			LRR-8	
	1,4	8	LTL-2	LTA-2	LTF-8		
		16	LRI-6			LRR-9	
	1,6	6		LRL-1		LRU-1	
	1,8	8	LTL-3	LTA-3	LTF-9	LRU-2	
	2,3	8	LTL-4	LRL-4	LTF-0	LRU-3	LRV-2 (6S)
		20	LRK-1	LRM-3		LRO-9	
	2,7	10	LTL-5	LTA-5	LTL-6		
	3,1	6	LRK-2	LTL-9	LRR-2	LRU-4	LRV-3
		10	LTK-1	LTA-6	LTF-1		
	3,5	6	LRK-4	LTL-7	LTF-2		
		12	LTK-2	LTA-7			
	4,0	6			LRG-8	LRR-3	LRV-4 (8S)
		12	LTK-3	LTA-8	LTF-3	LRU-6	
	4,5	6	LRK-7	LTL-8	LRO-1	LRR-4	
		12	LTK-4	LTA-9	LTF-4	LRU-7	
	5,0	8	LRK-9	LTL-0	LRO-2	LRU-8	LRV-5
		12	LTK-5	LTB-1	LTF-5		
	6,0	8		LRL-9		LRR-6	LRV-6
		14	LTK-6	LTB-2	LTF-6	LRU-9	
	7,0	8			LRO-4	LRR-7	
		16	LTK-7	LTB-3	LTF-7	LRV-1	


DIAMOND LOW SPEED Cutters / Fräser

						
	Ø mm	L 70 mm	L 95 mm	L 125 mm	L 55 mm	L 112 mm
	0,6	LUD-1	LUF-1	LUH-1		
	0,7	LUD-2		LUH-2		
	0,8	LUD-3	LUF-2	LUH-3		
	1,0	LUD-4	LUF-3	LUH-4	LRX-1	LRY-1
	1,2		LUF-4			
	1,4	LUD-5	LUF-5			
	1,6	LUD-0	LUF-6	LUH-5	LRX-2	
	1,8	LUD-6	LUF-7	LUH-6		
	2,1	LUD-7	LUF-8	LUH-7	LRX-3	
	2,3	LUD-8	LUF-9	LUH-8	LRX-4	LRY-2
	2,5	LUD-9		LUH-9		
	2,7	LUE-1	LUG-8	LUI-1		
	2,9		LUG-1		LRX-5	
	3,1	LUE-2	LUG-2	LUI-2		
	3,5	LUE-3	LUG-3	LUI-3	LRX-6	LRY-3
	4,0	LUE-4	LUG-4	LUI-4	LRX-7	LRY-4
	4,5	LUE-5	LUG-9	LUI-5	LRX-8	
	5,0	LUE-6	LUG-5	LUI-6	LRX-9	
	6,0	LUE-7	LUG-6	LUI-7		
	7,0	LUE-8	LUG-7	LUI-8		

CERAMO® LOW SPEED Cutters / Fräser

				
	Ø mm	L 70 mm	L 95 mm	L 125 mm
	4,0	LTH-4	LTC-1	LTD-1
	5,0	LTH-5	LTC-2	LTD-2
	6,0	LTH-6	LTC-3	LTD-3
	7,0	LTH-7	LTC-4	LTD-4
	8,0	LTH-8	LTC-5	LTD-5


CERAMO® LOW SPEED Twist Drills / Spiralbohrer



Zylindric Connection / Zylindrischer Anschluß

Ø mm	Length / Länge a mm	Length / Länge b mm	Item No. Art.-Nr.
1,1	35	20	LTY-1
1,5	60	40	LTY-2
2,0	75	50	LTY-3
2,5	135	50	LTY-5
2,5	75	50	LTY-6
2,5	95	50	LTY-7
2,7	75	50	LTY-8
3,2	145	50	LTY-4














CERAMO® LOW SPEED Twist Drills / Spiralbohrer











AO-Connection / AO-Kupplung

Ø mm	Length / Länge a mm	Length / Länge b mm	Item No. Art.-Nr.
1,1	60	20	LTN-5
1,5	85	30	LTN-6
1,5	110	30	LTO-2
1,5	125	30	LTO-9
1,8	100	30	LTN-7
2,0	100	50	LTN-8
2,0	125	50	LTP-1
2,4	100	50	LTN-9
2,5	180	50	LTO-7
2,5	110	50	LTO-3
2,7	100	50	LTO-1
2,7	125	50	LTP-2
3,2	195	50	LTP-3
3,2	145	50	LTO-5
3,2	225	50	LTP-7
3,5	110	50	LTO-4
3,5	195	50	LTP-4
4,0	225	50	LTP-8
4,5	145	50	LTO-6
4,5	195	50	LTP-5
6,0	195	50	LTP-6



















CERAMO® LINDEMANN Cutters / Fräser

	Ø mm	Length Business Part Länge Arbeitsteil	Total Length Gesamtlänge	Item No. Art.-Nr.
	1,4 mm	5,0	44,5	LRM-4
	1,4 mm	5,0	65,0	LRM-5
	1,4 mm	5,0	70,0	LRM-6
	1,6 mm	9,0	55,0	LRW-2
	2,3 mm	22,0	62,0	LRM-7
	2,3 mm	22,0	70,0	LRM-8
	2,3 mm	22,0	70,0	LRW-5
	2,3 mm	35,0	75,0	LRM-9
	2,3 mm	35,0	75,0	LRW-7
	2,3 mm	25,0	67,0	LRN-0
	2,3 mm	25,0	67,0	LRW-6
	2,3 mm	9,5	55,0	LRW-4













CERAMO® CRANIAL Cutters / Fräser

	L mm	Ø mm	equivalent to / äquivalent zu	Item No. Art.-Nr.
	44,5	2,0	3M C426 / Codman 26-1247	LRH-7
	44,5	2,5	Aesculap GC 303 / Codman 26-1246	LRH-8
	44,5	2,7	Z-Hall 5059-222	LRG-3
	44,5	2,7	3M C451	LRG-4
	44,5	2,9	Z-Hall 5059-220	LRG-5
	51,0	3,2	Z-Hall 5059-23	LRG-6
	51,0	2,9	3M C452	LRG-7
	51,0	3,7	Z-Hall 5059-25	LRG-0



















CERAMO® HIGH SPEED Rose Hobs / Rosenbohrer

									
	Ø mm	S	▲	I	II	III	O	M	AM
	1,0	6	LTJ-0	LTJ-7	LTM-4		LTX-3	LUB-0	LTT-3
	1,4	6	LTJ-1	LTJ-8	LTQ-0		LTX-4	LUB-1	LTT-4
	2,4	8	LTJ-2	LTJ-9	LTQ-1	LTQ-8	LTX-5	LUB-2	LTT-5
	3,2	10	LTJ-3	LTM-0	LTQ-2	LTQ-9	LTX-6	LUB-3	LTT-6
	4,0	8					LTQ-7	LUB-4	LTT-7
	4,8	10	LTJ-4	LTM-1	LTQ-3	LTX-0	LTX-7	LUB-5	LTT-8
	5,4	12						LUB-6	LTT-9
	6,4	14	LTJ-5	LTM-2	LTQ-4	LTX-1	LTX-8	LUB-7	LTS-5
	7,0	8						LUB-8	LTS-6
	7,9	16	LTJ-6	LTM-3	LTQ-5	LTX-2	LTX-9	LUB-9	LTS-7

CERAMO® HIGH SPEED Omni E Hobs / Omni E Bohrer

									
	Ø mm	S	▲	I	II	III	O	M	AM
	3,2	10		LUJ-0	LUJ-4			LUJ-8	LUK-2
	4,8	10		LUJ-1	LUJ-5			LUJ-9	LUK-3
	6,4	14		LUJ-2	LUJ-6			LUK-0	LUK-4
	7,9	16		LUJ-3	LUJ-7			LUK-1	LUK-5

DIAMOND HIGH SPEED Cutters / Fräser

								
	Ø mm	▲	I	II	III	O	M	AM
	1,0	LUK-6	LUU-4	LUV-1		LUW-1	LUW-8	LUX-8
	1,5	LUK-7	LUU-5	LUV-2		LUW-2	LUW-9	LUX-9
	2,1	LUK-8	LUU-6	LUV-3	LUV-7	LUW-3	LUX-0	LUY-0
	2,7					LUW-4	LUX-1	LUY-1
	3,3	LUK-9	LUU-7	LUV-4	LUV-8	LUW-5	LUX-2	LUY-2
	3,8		LUU-8				LUX-3	
	4,2					LUW-6	LUX-4	LUY-4
	4,6	LUU-2	LUU-9	LUV-5	LUV-9		LUX-5	
	5,0					LUW-7	LUX-6	LUY-6
	6,0	LUU-3	LUV-0	LUV-6	LUW-0		LUX-7	LUY-7