



CEBOTARI

UNIVERSAL STERNAL RETRACTOR UNIVERSAL STERNUM-SPREIZER




UNIVERSAL
STERNAL RETRACTOR
UNIVERSAL
STERNUM-SPREIZER

MQL-1



*exemplary configuration for IMA exposure
Konfigurationsbeispiel als IMA-Spreizer ¹*

*exemplary configuration with atrial retractor
Konfigurationsbeispiel mit atrialem Haken ²*



STERNOTOMY BLADE (PAIR)
STERNUMBLATT (PAAR)



IMA BLADE
IMA-BLATT

MQL-3

BLADE SIZE BLATTGRÖSSE a x b	
34 x 50 mm	MQL-4
43 x 50 mm	MQL-5
34 x 100 mm	MQL-6
43 x 100 mm	MQL-2
50 x 100 mm	MQL-7
63 x 100 mm	MQL-8
34 x 120 mm	MQL-9
43 x 120 mm	MQM-1
50 x 120 mm	MQM-2
63 x 120 mm	MQM-3

ACCESSORIES / ZUBEHÖR

¹ for IMA exposure / zur Darstellung der IMA

Baykut-IMA blade / Baykut-IMA-Blatt
MLC-2V
Catalogue / Katalog 2/8

² for atrial retraction / zur Vorhofretraktion

Clamping element / Befestigungselement +
Ball adapter / Kugeladapter +
Atrial retractor / Vorhofhaken
Catalogue / Katalog 2/16