



# FORMAR



PEEK CAGES		WIDTH BREITE	DEPTH TIEFE	HEIGHT HÖHE	TITANIUM SIZERS
STERILE	NON-STERILE				
PF-3S	PF-3	14 mm	13 mm	3 mm	SF-3
PF-4S	PF-4	14 mm	13 mm	4 mm	SF-4
PF-5S	PF-5	14 mm	13 mm	5 mm	SF-5
PF-6S	PF-6	14 mm	13 mm	6 mm	SF-6
PF-7S	PF-7	14 mm	13 mm	7 mm	SF-7
PG-3S	PG-3	16 mm	14 mm	3 mm	SG-3
PG-4S	PG-4	16 mm	14 mm	4 mm	SG-4
PG-5S	PG-5	16 mm	14 mm	5 mm	SG-5
PG-6S	PG-6	16 mm	14 mm	6 mm	SG-6
PG-7S	PG-7	16 mm	14 mm	7 mm	SG-7
PH-4S	PH-4	18 mm	14 mm	4 mm	SH-4
PH-5S	PH-5	18 mm	14 mm	5 mm	SH-5
PH-6S	PH-6	18 mm	14 mm	6 mm	SH-6
PH-7S	PH-7	18 mm	14 mm	7 mm	SH-7



DISSECTING CURETTE  
DISSEKTIONSKÜRETTE  
6 x 220 mm

XSE-1



STORAGE AND STERILISATION TRAY FOR CAGES + SIZERS  
LAGER- UND STERIBEHÄLTER FÜR CAGES + SIZERS



TITANIUM SETTING INSTRUMENT  
TITAN SETZ-INSTRUMENT  
180 mm

XRE-1

SMALL KLEIN	CAGE/SIZER TYPES CAGE/SIZER NUMMERN	PF + PG / SF + SG	PAZ-1
LARGE GROSS	CAGE/SIZER TYPES CAGE/SIZER NUMMERN	PF + PG + PH / SF + SG + SH	PAZ-2

### ALTERNATIVE TOOLS\* ALTERNATIVES WERKZEUG\*




TITANIUM SETTING INSTRUMENT WITH ADJUSTABLE STOP UNIT  
TITAN SETZ-INSTRUMENT MIT VERSTELLBAREM ANSCHLAG  
130/180 mm

XRE-4



TITANIUM SETTING INSTRUMENT WITH ADJUSTABLE STOP UNIT  
TITAN SETZ-INSTRUMENT MIT VERSTELLBAREM ANSCHLAG  
80/180 mm

XRE-3



CERAMO® LANKSCH  
INTERVERTEBRAL CURETTE  
INTERVERTEBRAL KÜRETTE  
L = 215 mm  
a = 100 mm

\*Does not fit into tray PAZ-1/-2

\*Passt nicht in Box PAZ-1/-2